

# BEAM SPLITTERS

## OP-TV

DOUBLE THE IMAGE IN THE EASIEST WAY



- Resistant
- Bright
- Focusable

OPTOMIC offers a complete range of video solutions for showing live images on a monitor, taking pictures or capturing videos, as well as the possibility of transferring the images to a PC or laptop.

The beam splitters are the perfect complement to our microscopes and colposcopes and are very easy to operate.

## BEAM SPLITTERS **OP-TV**



### ■ OP-TV CLICK

This beam splitter has a quick connection for easy attachment to any endoscopy camera through the endocoupler.

*A versatile solution for combining endoscopy and microscopy/colposcopy.*



### ■ OP-TV DC

Beam splitter for compact digital photographic camera (camera not included). Photos and videos with internal storage file. Consult compatible camera models.

*The best choice for taking still pictures and with the advantage of being wireless.*



### ■ OP-TV8

Compact system with integrated video camera.  
Video outputs: composite video and S-video (Y/C)

*An effective system for quality procedures.*



### ■ OP-TV10 HD

Compact system with integrated Full HD 1080p video camera.  
Video output: HDMI

*The professional solution with Full HD quality.*

**ALL OP-TV BEAM SPLITTERS ARE SIMPLE TO INSTALL AND EASY TO USE**

### **OPTOMIC ESPAÑA S.A.**

Madroño, 4. E-28770 Colmenar Viejo. Madrid, Spain

Tel: +34 91 804 44 44

[export@optomic.com](mailto:export@optomic.com)

[www.optomic.com](http://www.optomic.com)



producing quality

Specifications are subject to change without notice

